

产品名称: KI67 Mouse Monoclonal Antibody
产品货号: AMM80890



产品概述 (Summary)

产品名称 (Production Name)	KI67 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse,Rat,Rabbit,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KI67
别名 (Alternative Names)	KIA; Ki-67; MKI67
基因 ID (Gene ID)	4288.0
蛋白 ID (SwissProt ID)	P46013.Synthetic peptide corresponding to aa (CEDLAGFKELFQTPG) of human KI67, conjugated to KLH.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	358kDa

研究背景 (Background)

Ki67, also known as MKI67, it is the prototypic cell cycle related nuclear protein, expressed by proliferating cells in all


产品名称: KI67 Mouse Monoclonal Antibody
产品货号: AMM80890




phases of the active cell cycle (G1, S, G2 and M phase). It is absent in resting (G0) cells. Ki67 antibodies are useful in establishing the cell growing fraction in neoplasms (immunohistochemically quantified by determining the number of Ki67 positive cells among the total number of resting cells = Ki67 index). In neoplastic tissues the prognostic value is comparable to the tritiated thymidine labelling index. The correlation between low Ki67 index and histologically low grade tumours is strong. Ki67 is routinely used as a neuronal marker of cell cycling and proliferation.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis using KI67 mouse mAb against Hela (1), MCF-7 (2) and Raji (3) cell lysate.




Flow cytometric analysis of NIH/3T3 cells using KI67 mouse mAb (green) and negative control (red).


产品名称: KI67 Mouse Monoclonal Antibody
产品货号: AMM80890



Immunohistochemical analysis of paraffin-embedded human lung cancer (left) and rectal cancer (right) using KI67 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Mouse colon(left)Mouse (right) using KI67 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rabbit rectum(A) Rabbit spleen (B) Rabbit testis (C) Rabbit thymus (D) using KI67 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .