

产品名称: SNAI1 Mouse Monoclonal Antibody
产品货号: AMM80931



产品概述 (Summary)

产品名称 (Production Name)	SNAI1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SNAI1
别名 (Alternative Names)	SNA; SNAH; SLUGH2; dJ710H13.1; SNAI1
基因 ID (Gene ID)	6615.0
蛋白 ID (SwissProt ID)	O95863.Purified recombinant fragment of human SNAI1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	29kDa

研究背景 (Background)

Snail is a zinc-finger transcription factor that can repress E-cadherin transcription. Downregulation of E-cadherin is associated with epithelial-mesenchymal transition during embryonic development, a process also exploited by invasive


产品名称: SNAI1 Mouse Monoclonal Antibody
产品货号: AMM80931




cancer cells . Indeed, loss of E-cadherin expression is correlated with the invasive properties of some tumors and there is a considerable inverse correlation between Snail and E-cadherin mRNA levels in epithelial tumor cell lines . In addition, Snail blocks the cell cycle and confers resistance to cell death . Phosphorylation of Snail by GSK-3 and PAK1 regulates its stability, cellular localization and function .Tissue specificity: Expressed in a variety of tissues with the highest expression in kidney.

研究领域 (Research Area)

图片 (Image Data)



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Western blot analysis using SNAI1 mouse mAb against NTERA-2 cell lysate.

注意事项 (Note)

For research use only .