

产品名称: RBP4 Mouse Monoclonal Antibody
产品货号: AMM80976



产品概述 (Summary)

产品名称 (Production Name)	RBP4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBP4
别名 (Alternative Names)	RBP4
基因 ID (Gene ID)	5950.0
蛋白 ID (SwissProt ID)	P02753. Purified recombinant fragment of human RBP expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:200-1:1000, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	23kDa

研究背景 (Background)

This protein belongs to the lipocalin family and is the specific carrier for retinol (vitamin A alcohol) in the blood. It delivers retinol from the liver stores to the peripheral tissues. In plasma, the RBP-retinol complex interacts with transthyretin which


产品名称: RBP4 Mouse Monoclonal Antibody
产品货号: AMM80976




prevents its loss by filtration through the kidney glomeruli. A deficiency of vitamin A blocks secretion of the binding protein posttranslationally and results in defective delivery and supply to the epidermal cells. (provided by RefSeq)

研究领域 (Research Area)


图片 (Image Data)



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);




Immunofluorescence analysis of HepG2 cells using RBP4 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow cytometric analysis of HepG2 cells using RBP4 mouse mAb (green) and negative control (purple).

产品名称: RBP4 Mouse Monoclonal Antibody
产品货号: AMM80976



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues (left) and stomach cancer tissues (right) using RBP4 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .