

产品名称: Epha10 Mouse Monoclonal Antibody
产品货号: AMM80984



产品概述 (Summary)

产品名称 (Production Name)	Epha10 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Epha10
别名 (Alternative Names)	FLJ16103; FLJ33655; MGC43817
基因 ID (Gene ID)	284656.0
蛋白 ID (SwissProt ID)	Q5JZY3. Purified recombinant fragment of human Epha10 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	110kDa

研究背景 (Background)

Ephrin receptors, the largest subfamily of receptor tyrosine kinases (RTKs), and their ephrin ligands are important mediators of cell-cell communication regulating cell attachment, shape, and mobility in neuronal and epithelial cells (Aasheim et al.,


产品名称: Epha10 Mouse Monoclonal Antibody
产品货号: AMM80984



2005 (PubMed 15777695)).Tissue specificity: Mainly expressed in testis.

研究领域 (Research Area)

图片 (Image Data)



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);

注意事项 (Note)

For research use only .