

产品名称: ApoE Mouse Monoclonal Antibody
产品货号: AMM80990



产品概述 (Summary)

产品名称 (Production Name)	ApoE Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ApoE
别名 (Alternative Names)	AD2; LPG; LDLCQ5; MGC1571
基因 ID (Gene ID)	348.0
蛋白 ID (SwissProt ID)	P02649.Purified recombinant fragment of human ApoE expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	36kDa

研究背景 (Background)

Chylomicron remnants and very low density lipoprotein (VLDL) remnants are rapidly removed from the circulation by receptor-mediated endocytosis in the liver. Apolipoprotein E, a main apoprotein of the chylomicron, binds to a specific


产品名称: ApoE Mouse Monoclonal Antibody
产品货号: AMM80990




receptor on liver cells and peripheral cells. ApoE is essential for the normal catabolism of triglyceride-rich lipoprotein constituents. The APOE gene is mapped to chromosome 19 in a cluster with APOC1 and APOC2. Defects in apolipoprotein E result in familial dysbetalipoproteinemia, or type III hyperlipoproteinemia (HLP III), in which increased plasma cholesterol and triglycerides are the consequence of impaired clearance of chylomicron and VLDL remnants. Tissue specificity: Occurs in all lipoprotein fractions in plasma. It constitutes 10-20% of very low density lipoproteins (VLDL) and 1-2% of high density lipoproteins (HDL). APOE is produced in most organs. Significant quantities are produced in liver, brain, spleen, lung, adrenal, ovary, kidney and muscle.

研究领域 (Research Area)

图片 (Image Data)




Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Flow cytometric analysis of HepG2 cells using ApoE mouse mAb (green) and negative control (purple).

产品名称: ApoE Mouse Monoclonal Antibody
产品货号: AMM80990



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues (left) and brain tissues (right) using ApoE mouse mAb with DAB staining.

注意事项 (Note)

For research use only .