

产品名称: NFKB1 Mouse Monoclonal Antibody
产品货号: AMM80999



产品概述 (Summary)

产品名称 (Production Name)	NFKB1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NFKB1
别名 (Alternative Names)	p50; KBF1; p105; EBP-1; MGC54151; NFKB-p50; NFkappaB; NF-kappaB; NFKB-p105; NF-kappa-B
基因 ID (Gene ID)	4790.0
蛋白 ID (SwissProt ID)	P19838.Purified recombinant fragment of human NFKB1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	50kDa/105kDa

研究背景 (Background)

This gene encodes a 105 kD protein which can undergo cotranslational processing by the 26S proteasome to produce a 50

产品名称: NFKB1 Mouse Monoclonal Antibody
产品货号: AMM80999




kD protein. The 105 kD protein is a Rel protein-specific transcription inhibitor and the 50 kD protein is a DNA binding subunit of the NF-kappa-B (NFKB) protein complex. NFKB is a transcription regulator that is activated by various intra- and extra-cellular stimuli such as cytokines, oxidant-free radicals, ultraviolet irradiation, and bacterial or viral products. Activated NFKB translocates into the nucleus and stimulates the expression of genes involved in a wide variety of biological functions. Inappropriate activation of NFKB has been associated with a number of inflammatory diseases while persistent inhibition of NFKB leads to inappropriate immune cell development or delayed cell growth. Two transcript variants encoding different isoforms have been found for this gene.


研究领域 (Research Area)

Apoptosis, PI3K-Akt signaling pathway, MAPK signaling pathway

图片 (Image Data)




Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);




Western blot analysis using NFKB1 mouse mAb against K562 (1), Jurkat (2), A431 (3), HeLa (4), THP-1 (5) and MCF-7 (6) cell lysate.

产品名称: NFKB1 Mouse Monoclonal Antibody
产品货号: AMM80999



Flow cytometric analysis of MCF-7 cells using NFKB1 mouse mAb (green) and negative control (purple).



Immunohistochemical analysis of paraffin-embedded human prostate tissues (left) and bladder cancer tissues (right) using NFKB1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .