

产品名称: CRKL Mouse Monoclonal Antibody
产品货号: AMM81006



产品概述 (Summary)

产品名称 (Production Name)	CRKL Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CRKL
别名 (Alternative Names)	CRKL
基因 ID (Gene ID)	1399.0
蛋白 ID (SwissProt ID)	P46109.Purified recombinant fragment of human CRKL expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	39kDa

研究背景 (Background)

This gene encodes a protein kinase containing SH2 and SH3 (src homology) domains which has been shown to activate the RAS and JUN kinase signaling pathways and transform fibroblasts in a RAS-dependent fashion. It is a substrate of the BCR-

产品名称: CRKL Mouse Monoclonal Antibody
产品货号: AMM81006




ABL tyrosine kinase, plays a role in fibroblast transformation by BCR-ABL, and may be oncogenic.


研究领域 (Research Area)

MAPK signaling pathway

图片 (Image Data)



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Flow cytometric analysis of NIH/3T3 cells using CRKL mouse mAb (blue) and negative control (red).

注意事项 (Note)

For research use only .