

产品名称: GKAP Mouse Monoclonal Antibody
产品货号: AMM81013



产品概述 (Summary)

产品名称 (Production Name)	GKAP Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ICC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GKAP
别名 (Alternative Names)	GKAP; DLGAP1; DAP-1; hGKAP; SAPAP1; FLJ38442; MGC88156; DAP-1-BETA; DAP-1-ALPHA
基因 ID (Gene ID)	9229.0
蛋白 ID (SwissProt ID)	O14490.Purified recombinant fragment of human GKAP expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	109kDa

研究背景 (Background)


Function: Part of the postsynaptic scaffold in neuronal cells.Tissue specificity: Expressed in brain.

产品名称: GKAP Mouse Monoclonal Antibody
产品货号: AMM81013




研究领域 (Research Area)


图片 (Image Data)



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);




Immunofluorescence analysis of NIH/3T3 cells using GKAP mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using GKAP mouse mAb with DAB staining.

产品名称: GKAP Mouse Monoclonal Antibody
产品货号: AMM81013



Immunohistochemical analysis of paraffin-embedded human esophagus tissues using GKAP mouse mAb with DAB staining.

注意事项 (Note)

For research use only .