

产品名称: CRK Mouse Monoclonal Antibody
产品货号: AMM81030



产品概述 (Summary)

产品名称 (Production Name)	CRK Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CRK
别名 (Alternative Names)	CRKII
基因 ID (Gene ID)	1398.0
蛋白 ID (SwissProt ID)	P46108. Purified recombinant fragment of human CRK expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:200-1:1000, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	42kDa

研究背景 (Background)

This gene encodes a member of an adapter protein family that binds to several tyrosine-phosphorylated proteins. The product of this gene has several SH2 and SH3 domains (src-homology domains) and is involved in several signaling

产品名称: CRK Mouse Monoclonal Antibody
产品货号: AMM81030




pathways, recruiting cytoplasmic proteins in the vicinity of tyrosine kinase through SH2-phosphotyrosine interaction. The N-terminal SH2 domain of this protein functions as a positive regulator of transformation whereas the C-terminal SH3 domain functions as a negative regulator of transformation. Two alternative transcripts encoding different isoforms with distinct biological activity have been described.


研究领域 (Research Area)

MAPK signaling pathway

图片 (Image Data)




Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);




Immunofluorescence analysis of 3T3-L1 cells using CRK mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.


产品名称: CRK Mouse Monoclonal Antibody
产品货号: AMM81030



Flow cytometric analysis of Hela cells using CRK mouse mAb (blue) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human intima cancer tissues using CRK mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using CRK mouse mAb with DAB staining.

注意事项 (Note)

For research use only .