

产品名称: HFE Mouse Monoclonal Antibody
产品货号: AMM81041



产品概述 (Summary)

产品名称 (Production Name)	HFE Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HFE
别名 (Alternative Names)	HH; HFE1; HLA-H; MVCD7; TFQTL2; MGC103790; dJ221C16.10.1
基因 ID (Gene ID)	3077.0
蛋白 ID (SwissProt ID)	Q30201.Purified recombinant fragment of human HFE expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	40kDa

研究背景 (Background)

The protein encoded by this gene is a membrane protein that is similar to MHC class I-type proteins and associates with beta2-microglobulin (beta2M). It is thought that this protein functions to regulate iron absorption by regulating the


产品名称: HFE Mouse Monoclonal Antibody
产品货号: AMM81041




interaction of the transferrin receptor with transferrin. The iron storage disorder, hereditary haemochromatosis, is a recessive genetic disorder that results from defects in this gene. At least nine alternatively spliced variants have been described for this gene. Additional variants have been found but their full-length nature has not been determined. bsp;

研究领域 (Research Area)

图片 (Image Data)



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Immunofluorescence analysis of HepG2 cells using HFE mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.

注意事项 (Note)

For research use only .