

产品名称: TNNI2 Mouse Monoclonal Antibody
产品货号: AMM81045



产品概述 (Summary)

产品名称 (Production Name)	TNNI2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNNI2
别名 (Alternative Names)	DA2B; FSSV; fsTnl; AMCD2B
基因 ID (Gene ID)	7136.0
蛋白 ID (SwissProt ID)	P48788.Purified recombinant fragment of human TNNI2 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	21kDa

研究背景 (Background)

This gene encodes a fast-twitch skeletal muscle protein, a member of the troponin I gene family, and a component of the troponin complex including troponin T, troponin C and troponin I subunits. The troponin complex, along with tropomyosin,


产品名称: TNNI2 Mouse Monoclonal Antibody
产品货号: AMM81045




is responsible for the calcium-dependent regulation of striated muscle contraction. Mouse studies show that this component is also present in vascular smooth muscle and may play a role in regulation of smooth muscle function. In addition to muscle tissues, this protein is found in corneal epithelium, cartilage where it is an inhibitor of angiogenesis to inhibit tumor growth and metastasis, and mammary gland where it functions as a co-activator of estrogen receptor-related receptor alpha. This protein also suppresses tumor growth in human ovarian carcinoma. Mutations in this gene cause myopathy and distal arthrogryposis type 2B. Alternatively spliced transcript variants have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)




Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);




Flow cytometric analysis of NIH/3T3 cells using TNNI2 mouse mAb (blue) and negative control (red).


产品名称: TNNI2 Mouse Monoclonal Antibody
产品货号: AMM81045



Immunohistochemical analysis of paraffin-embedded human rabbit cardiac muscle tissues using TNNI2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human cardiac muscle tissues using TNNI2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human striated muscle tissues using TNNI2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .