

产品名称: LPlunc1 Mouse Monoclonal Antibody
产品货号: AMM81071



产品概述 (Summary)

产品名称 (Production Name)	LPlunc1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LPlunc1
别名 (Alternative Names)	BPIFB1; LPLUNC1; MGC14597; C20orf114
基因 ID (Gene ID)	92747.0
蛋白 ID (SwissProt ID)	Q8TDL5. Purified recombinant fragment of human LPlunc1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:200-1:1000, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	52kDa

研究背景 (Background)

The protein encoded by this gene may be involved in the innate immune response to bacterial exposure in the mouth, nasal cavities, and lungs. The encoded protein is secreted and is a member of the BPI/LBP/PLUNC protein superfamily. This


产品名称: LPlunc1 Mouse Monoclonal Antibody
产品货号: AMM81071




gene is found with other members of the superfamily in a cluster on chromosome 20.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Immunofluorescence analysis of HeLa cells using LPlunc1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.




Flow cytometric analysis of HeLa cells using LPlunc1 mouse mAb (green) and negative control (red).

产品名称: LPlunc1 Mouse Monoclonal Antibody
产品货号: AMM81071



Immunohistochemical analysis of paraffin-embedded human brain tissues using LPlunc1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human tonsil tissues using LPlunc1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .