

产品名称: THPO Mouse Monoclonal Antibody
产品货号: AMM81119



产品概述 (Summary)

产品名称 (Production Name)	THPO Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	THPO
别名 (Alternative Names)	ML; TPO; MGDF; MKCSF; MPLLG
基因 ID (Gene ID)	7066.0
蛋白 ID (SwissProt ID)	P40225.Purified recombinant fragment of human THPO expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	37.8kDa

研究背景 (Background)

Megakaryocytopoiesis is the cellular development process that leads to platelet production. The protein encoded by this gene is a humoral growth factor that is necessary for megakaryocyte proliferation and maturation, as well as for

产品名称: THPO Mouse Monoclonal Antibody
产品货号: AMM81119




thrombopoiesis. This protein is the ligand for MLP/C_MPL, the product of myeloproliferative leukemia virus oncogene. Alternate splicing results in multiple transcript variants of this gene.


研究领域 (Research Area)

TGF-beta signaling pathway

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of MCF-7 cells using THPO mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .