

产品名称: HAS3 Mouse Monoclonal Antibody
产品货号: AMM81124



产品概述 (Summary)

产品名称 (Production Name)	HAS3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HAS3
别名 (Alternative Names)	HAS3
基因 ID (Gene ID)	3038.0
蛋白 ID (SwissProt ID)	O00219 .Purified recombinant fragment of human HAS3 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	63kDa

研究背景 (Background)

The protein encoded by this gene is involved in the synthesis of the unbranched glycosaminoglycan hyaluronan, or hyaluronic acid, which is a major constituent of the extracellular matrix. This gene is a member of the NODC/HAS gene


产品名称: HAS3 Mouse Monoclonal Antibody
产品货号: AMM81124




family. Compared to the proteins encoded by other members of this gene family, this protein appears to be more of a regulator of hyaluronan synthesis. Alternative splicing results in multiple transcript variants.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Flow cytometric analysis of HeLa cells using HAS3 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using HAS3 mouse mAb with DAB staining.

产品名称: HAS3 Mouse Monoclonal Antibody
产品货号: AMM81124



Immunohistochemical analysis of paraffin-embedded human prostate tissues using HAS3 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .