

产品名称: GAB1 Mouse Monoclonal Antibody
产品货号: AMM81125



产品概述 (Summary)

产品名称 (Production Name)	GAB1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GAB1
别名 (Alternative Names)	GAB1
基因 ID (Gene ID)	2549.0
蛋白 ID (SwissProt ID)	Q13480.Purified recombinant fragment of human GAB1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	76.6kDa

研究背景 (Background)

The protein encoded by this gene is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis.


产品名称: GAB1 Mouse Monoclonal Antibody
产品货号: AMM81125




Two transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using GAB1 mouse mAb against HEK293 (1) cell lysate.




Flow cytometric analysis of Jurkat cells using GAB1 mouse mAb (green) and negative control (red).

产品名称: GAB1 Mouse Monoclonal Antibody
产品货号: AMM81125



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using GAB1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human kidney tissues using GAB1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .