

产品名称: mouse Splunc2 Mouse Monoclonal Antibody
产品货号: AMM81126



产品概述 (Summary)

产品名称 (Production Name)	mouse Splunc2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Splunc2
别名 (Alternative Names)	Splunc2
基因 ID (Gene ID)	19194.0
蛋白 ID (SwissProt ID)	P07743 .Purified recombinant fragment of mouse Splunc2 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	/

研究背景 (Background)


N/A

产品名称: mouse Splunc2 Mouse Monoclonal Antibody
产品货号: AMM81126




研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Immunofluorescence analysis of HepG2 cells. Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.




Immunofluorescence analysis of HepG2 cells using mouse Splunc2 mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.


产品名称: mouse Splunc2 Mouse Monoclonal Antibody
产品货号: AMM81126



Flow cytometric analysis of HeLa cells using mouse Splunc2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human prostate tissues using mouse Splunc2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human kidney tissues using mouse Splunc2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .