

产品名称: NCK1 Mouse Monoclonal Antibody
产品货号: AMM81129



产品概述 (Summary)

产品名称 (Production Name)	NCK1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCK1
别名 (Alternative Names)	NCK; nck-1; NCKalpha
基因 ID (Gene ID)	4690.0
蛋白 ID (SwissProt ID)	P16333.Purified recombinant fragment of human NCK1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	47kDa

研究背景 (Background)

The protein encoded by this gene is one of the signaling and transforming proteins containing Src homology 2 and 3 (SH2 and SH3) domains. It is located in the cytoplasm and is an adaptor protein involved in transducing signals from receptor


产品名称: NCK1 Mouse Monoclonal Antibody
产品货号: AMM81129




tyrosine kinases to downstream signal recipients such as RAS. Alternatively spliced transcript variants encoding different isoforms have been found.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using NCK1 mouse mAb against Jurkat (1), HeLa (2), HEK293 (3), A431 (4), K562 (5), and COS7 (6) cell lysate.



Flow cytometric analysis of Jurkat cells using NCK1 mouse mAb (green) and negative control (red).

产品名称: NCK1 Mouse Monoclonal Antibody
产品货号: AMM81129



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using NCK1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .