

产品名称: DLL4 Mouse Monoclonal Antibody
产品货号: AMM81140



产品概述 (Summary)

产品名称 (Production Name)	DLL4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DLL4
别名 (Alternative Names)	hdelta2
基因 ID (Gene ID)	54567.0
蛋白 ID (SwissProt ID)	Q9NR61.Purified recombinant fragment of human DLL4 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	74.6kDa

研究背景 (Background)

This gene is a homolog of the Drosophila delta gene. The delta gene family encodes Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain.


产品名称: DLL4 Mouse Monoclonal Antibody
产品货号: AMM81140




研究领域 (Research Area)

Notch signaling pathway

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HeLa cells using DLL4 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .