

产品名称: PPM1A Mouse Monoclonal Antibody
产品货号: AMM81143



产品概述 (Summary)

产品名称 (Production Name)	PPM1A Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPM1A
别名 (Alternative Names)	PP2CA; PP2Calpha; PP2C-ALPHA
基因 ID (Gene ID)	5494.0
蛋白 ID (SwissProt ID)	P20265 .Purified recombinant fragment of human PPM1A (AA: 202-382) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	42.4kDa

研究背景 (Background)

The protein encoded by this gene is a member of the PP2C family of Ser/Thr protein phosphatases. PP2C family members


产品名称: PPM1A Mouse Monoclonal Antibody
产品货号: AMM81143




are known to be negative regulators of cell stress response pathways. This phosphatase dephosphorylates, and negatively regulates the activities of, MAP kinases and MAP kinase kinases. It has been shown to inhibit the activation of p38 and JNK kinase cascades induced by environmental stresses. This phosphatase can also dephosphorylate cyclin-dependent kinases, and thus may be involved in cell cycle control. Overexpression of this phosphatase is reported to activate the expression of the tumor suppressor gene TP53/p53, which leads to G2/M cell cycle arrest and apoptosis. Three alternatively spliced transcript variants encoding distinct isoforms have been described.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using PPM1A mouse mAb against Jurkat (1), Jurkat (2), A431 (3), HeLa (4), HEK293 (5), Raji (6), MCF-7 (7), and COS7 (8) cell lysate.

产品名称: PPM1A Mouse Monoclonal Antibody
产品货号: AMM81143



Flow cytometric analysis of HeLa cells using PPM1A mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .