

产品名称: MAP2 Mouse Monoclonal Antibody
产品货号: AMM81144



产品概述 (Summary)

产品名称 (Production Name)	MAP2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAP2
别名 (Alternative Names)	MAP2A; MAP2B; MAP2C
基因 ID (Gene ID)	4133.0
蛋白 ID (SwissProt ID)	P11137.Purified recombinant fragment of human MAP2 (AA: 24-123) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	199kDa

研究背景 (Background)

This gene encodes a protein that belongs to the microtubule-associated protein family. The proteins of this family are


产品名称: MAP2 Mouse Monoclonal Antibody
产品货号: AMM81144




thought to be involved in microtubule assembly, which is an essential step in neurogenesis. The products of similar genes in rat and mouse are neuron-specific cytoskeletal proteins that are enriched in dendrites, implicating a role in determining and stabilizing dendritic shape during neuron development. A number of alternatively spliced variants encoding distinct isoforms have been described.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HeLa cells using MAP2 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .