

产品名称: PAX5 Mouse Monoclonal Antibody
产品货号: AMM81161



产品概述 (Summary)

产品名称 (Production Name)	PAX5 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PAX5
别名 (Alternative Names)	BSAP
基因 ID (Gene ID)	5079.0
蛋白 ID (SwissProt ID)	Q02548.Purified recombinant fragment of human PAX5 (AA: 235-382) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	45kDa

研究背景 (Background)

This gene encodes a member of the paired box (PAX) family of transcription factors. The central feature of this gene family


产品名称: PAX5 Mouse Monoclonal Antibody
产品货号: AMM81161




is a novel, highly conserved DNA-binding motif, known as the paired box. PAX proteins are important regulators in early development, and alterations in the expression of their genes are thought to contribute to neoplastic transformation. This gene encodes the B-cell lineage specific activator protein that is expressed at early, but not late stages of B-cell differentiation. Its expression has also been detected in developing CNS and testis and so the encoded protein may also play a role in neural development and spermatogenesis. This gene is located at 9p13, which is involved in t(9;14)(p13;q32) translocations recurring in small lymphocytic lymphomas of the plasmacytoid subtype, and in derived large-cell lymphomas. This translocation brings the potent E-mu enhancer of the IgH gene into close proximity of the PAX5 promoter, suggesting that the deregulation of transcription of this gene contributes to the pathogenesis of these lymphomas. Alternatively spliced transcript variants encoding different isoforms have been described but their biological validity has not been determined.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human brain tumors tissues using PAX5 mouse mAb with DAB staining.

产品名称: PAX5 Mouse Monoclonal Antibody
产品货号: AMM81161



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using PAX5 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .