

产品名称: SPIB Mouse Monoclonal Antibody  
产品货号: AMM81169



## 产品概述 (Summary)

|                        |                                |
|------------------------|--------------------------------|
| 产品名称 (Production Name) | SPIB Mouse Monoclonal Antibody |
| 描述 (Description)       | Mouse monoclonal Antibody      |
| 宿主 (Host)              | Mouse                          |
| 应用 (Application)       | WB,IHC,ICC,ELISA,FC            |
| 种属反应性 (Reactivity)     | Human,Mouse                    |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | Mouse IgG1                                                                               |
| 克隆 (Clonality)      | Monoclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide.                                        |
| 纯化方式 (Purification) | Affinity Purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                                                         |
|------------------------|-----------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | SPIB                                                                                    |
| 别名 (Alternative Names) | SPI-B                                                                                   |
| 基因 ID (Gene ID)        | 6689.0                                                                                  |
| 蛋白 ID (SwissProt ID)   | Q01892 .Purified recombinant fragment of human SPIB (AA: 200-252) expressed in E. Coli. |

## 产品应用 (Application)

|                          |                                                                                       |
|--------------------------|---------------------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400 |
| 蛋白分子量 (Molecular Weight) | 28.8kDa                                                                               |

## 研究背景 (Background)


产品名称: SPIB Mouse Monoclonal Antibody  
产品货号: AMM81169




The protein encoded by this gene is a transcriptional activator that binds to the PU-box (5'-GAGGAA-3') and acts as a lymphoid-specific enhancer. Four transcript variants encoding different isoforms have been found for this gene.

## 研究领域 (Research Area)


## 图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using SPIB mouse mAb against A549 (1), PC-3 (2), and NIH3T3 (3) cell lysate.




Immunofluorescence analysis of HeLa cells using SPIB mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.


产品名称: SPIB Mouse Monoclonal Antibody  
产品货号: AMM81169



Flow cytometric analysis of NIH3T3 cells using SPIB mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using SPIB mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer with DAB staining.

### 注意事项 (Note)

For research use only .