

产品名称: FZD5 Mouse Monoclonal Antibody
产品货号: AMM81177



产品概述 (Summary)

产品名称 (Production Name)	FZD5 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FZD5
别名 (Alternative Names)	HFZ5; C2orf31; MGC129692; DKFZp434E2135
基因 ID (Gene ID)	7855.0
蛋白 ID (SwissProt ID)	Q13467.Purified recombinant fragment of human FZD5 (AA:151-217) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	64.5kDa

研究背景 (Background)

Members of the 'frizzled' gene family encode 7-transmembrane domain proteins that are receptors for Wnt signaling

产品名称: FZD5 Mouse Monoclonal Antibody
产品货号: AMM81177




proteins. The FZD5 protein is believed to be the receptor for the Wnt5A ligand.


研究领域 (Research Area)

Wnt signaling pathway


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using FZD5 mouse mAb against A549 (1), and PC-3 (2) cell lysate.



Flow cytometric analysis of HepG2 cells using FZD5 mouse mAb (green) and negative control (red).

产品名称: FZD5 Mouse Monoclonal Antibody
产品货号: AMM81177



注意事项 (Note)

For research use only .