

产品名称: SCGB1A1 Mouse Monoclonal Antibody
产品货号: AMM81190



产品概述 (Summary)

产品名称 (Production Name)	SCGB1A1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SCGB1A1
别名 (Alternative Names)	UGB; UP1; CC10; CC16; CCSP; UP-1; CCPBP
基因 ID (Gene ID)	7356.0
蛋白 ID (SwissProt ID)	P11684.Purified recombinant fragment of human SCGB1A1 (AA: 26-91) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	10kDa

研究背景 (Background)

This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been


产品名称: SCGB1A1 Mouse Monoclonal Antibody
产品货号: AMM81190




implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of A549 cells using SCGB1A1 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .