

产品名称: KEAP1 Mouse Monoclonal Antibody
产品货号: AMM81199



产品概述 (Summary)

产品名称 (Production Name)	KEAP1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse,Rat,Rabbit,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KEAP1
别名 (Alternative Names)	INrf2; KLHL19
基因 ID (Gene ID)	9817.0
蛋白 ID (SwissProt ID)	Q14145.Purified recombinant fragment of human KEAP1 (AA: 380-624) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	69.7kDa

研究背景 (Background)

This gene encodes a protein containing KELCH-1 like domains, as well as a BTB/POZ domain. Kelch-like ECH-associated


产品名称: KEAP1 Mouse Monoclonal Antibody
产品货号: AMM81199




protein 1 interacts with NF-E2-related factor 2 in a redox-sensitive manner and the dissociation of the proteins in the cytoplasm is followed by transportation of NF-E2-related factor 2 to the nucleus. This interaction results in the expression of the catalytic subunit of gamma-glutamylcysteine synthetase. Two alternatively spliced transcript variants encoding the same isoform have been found for this gene.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Flow cytometric analysis of HeLa cells using KEAP1 mouse mAb (green) and negative control (red).




产品名称: KEAP1 Mouse Monoclonal Antibody
产品货号: AMM81199




Flow cytometric analysis of NIH/3T3 cells using KEAP1 mouse mAb (green) and negative control (red).



Flow cytometric analysis of COS7 cells using KEAP1 mouse mAb (green) and negative control (red).




Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using KEAP1 mouse mAb with DAB staining.




Immunohistochemical analysis of paraffin-embedded human cerebellum tissues using KEAP1 mouse mAb with DAB staining.


产品名称: KEAP1 Mouse Monoclonal Antibody
产品货号: AMM81199



Immunohistochemical analysis of paraffin-embedded Mouse brain using KEAP1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rat brain using KEAP1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rat muscle using KEAP1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .