

产品名称: WDR66 Mouse Monoclonal Antibody
产品货号: AMM81209



产品概述 (Summary)

产品名称 (Production Name)	WDR66 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	WDR66
别名 (Alternative Names)	WDR66
基因 ID (Gene ID)	144406.0
蛋白 ID (SwissProt ID)	Q8TBY9.Purified recombinant fragment of human WDR66 (AA: 1-250) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	130kDa

研究背景 (Background)

This protein encoded by this gene belongs to the WD repeat-containing family of proteins, which function in the formation


产品名称: WDR66 Mouse Monoclonal Antibody
产品货号: AMM81209




of protein-protein complexes in a variety of biological pathways. This family member appears to function in the determination of mean platelet volume (MPV), and polymorphisms in this gene have been associated with variance in MPV. Alternative splicing of this gene results in multiple transcript variants.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HepG2 cells using WDR66 mouse mAb (green) and negative control (purple).

注意事项 (Note)

For research use only .