

产品名称: GATA6 Mouse Monoclonal Antibody
产品货号: AMM81213



产品概述 (Summary)

产品名称 (Production Name)	GATA6 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GATA6
别名 (Alternative Names)	GATA6
基因 ID (Gene ID)	2627.0
蛋白 ID (SwissProt ID)	Q92908.Purified recombinant fragment of human GATA6 (AA: 491-557) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	60kDa

研究背景 (Background)

This gene is a member of a small family of zinc finger transcription factors that play an important role in the regulation of


产品名称: GATA6 Mouse Monoclonal Antibody
产品货号: AMM81213




cellular differentiation and organogenesis during vertebrate development. This gene is expressed during early embryogenesis and localizes to endo- and mesodermally derived cells during later embryogenesis and thereby plays an important role in gut, lung, and heart development. Mutations in this gene are associated with several congenital defects.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HEK293 cells using GATA6 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .