

产品名称: GSC Mouse Monoclonal Antibody
产品货号: AMM81223



产品概述 (Summary)

产品名称 (Production Name)	GSC Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSC
别名 (Alternative Names)	GSC
基因 ID (Gene ID)	145258.0
蛋白 ID (SwissProt ID)	P56915 .Purified recombinant fragment of human GSC (AA: 191-257) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	28.2kDa

研究背景 (Background)

This gene encodes a member of the bicoid subfamily of the paired (PRD) homeobox family of proteins. The encoded


产品名称: GSC Mouse Monoclonal Antibody
产品货号: AMM81223




protein acts as a transcription factor and may be autoregulatory. A similar protein in mice plays a role in craniofacial and rib cage development during embryogenesis.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Flow cytometric analysis of Hela cells using GSC mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using GSC mouse mAb with DAB staining.

产品名称: GSC Mouse Monoclonal Antibody
产品货号: AMM81223



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using GSC mouse mAb with DAB staining.

注意事项 (Note)

For research use only .