

产品名称: PLAGL1 Mouse Monoclonal Antibody
产品货号: AMM81246



产品概述 (Summary)

产品名称 (Production Name)	PLAGL1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PLAGL1
别名 (Alternative Names)	ZAC; LOT1; ZAC1
基因 ID (Gene ID)	5325.0
蛋白 ID (SwissProt ID)	Q9UM63.Purified recombinant fragment of human PLAGL1 (AA: 118-222) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	50.8kDa

研究背景 (Background)

This gene encodes a C2H2 zinc finger protein with transactivation and DNA-binding activities. It has been shown to have


产品名称: PLAGL1 Mouse Monoclonal Antibody
产品货号: AMM81246




anti-proliferative properties, and thus thought to function as a tumor suppressor. In addition, overexpression of this gene during fetal development is believed to underlie the rare disorder, transient neonatal diabetes mellitus (TNDM). This gene is imprinted, with preferential expression of the paternal allele in many tissues, however, biallelic expression has been noted in peripheral blood leucocytes. A recent study reports that tissue-specific imprinting results from variable utilization of monoallelic and biallelic promoters. Many transcript variants differing in the 5' UTR and encoding two different isoforms, have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using PLAGL1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .