

产品名称: ONECUT3 Mouse Monoclonal Antibody
产品货号: AMM81249



产品概述 (Summary)

产品名称 (Production Name)	ONECUT3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ICC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ONECUT3
别名 (Alternative Names)	OC3
基因 ID (Gene ID)	390874.0
蛋白 ID (SwissProt ID)	O60422.Purified recombinant fragment of human ONECUT3 (AA: 350-455) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	50kDa

研究背景 (Background)

The protein encoded by this gene is a transcriptional activator. It can bind the consensus DNA sequence 5'-


产品名称: ONECUT3 Mouse Monoclonal Antibody
产品货号: AMM81249




DHWATTGAYTWW-3' on a variety of gene promoters such as those of HNF3B and TTR(By similarity).
/>

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Immunofluorescence analysis of A431 cells using ONECUT3 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using ONECUT3 mouse mAb with DAB staining.

产品名称: ONECUT3 Mouse Monoclonal Antibody
产品货号: AMM81249



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using ONECUT3 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .