

产品名称: DNMT1 Mouse Monoclonal Antibody
产品货号: AMM81255



产品概述 (Summary)

产品名称 (Production Name)	DNMT1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNMT1
别名 (Alternative Names)	AIM; DNMT; MCMT; CXXC9; HSN1E
基因 ID (Gene ID)	1786.0
蛋白 ID (SwissProt ID)	P26358.Purified recombinant fragment of human DNMT1 (AA: 1448-1594) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	183.2kDa

研究背景 (Background)

DNA (cytosine-5-)-methyltransferase 1 has a role in the establishment and regulation of tissue-specific patterns of


产品名称: DNMT1 Mouse Monoclonal Antibody
产品货号: AMM81255




methyated cytosine residues. Aberrant methylation patterns are associated with certain human tumors and developmental abnormalities. Two transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using IL2RA mouse mAb against Jurkat (1), Cos7 (2), HCT116 (3) and NTERA-2 (4) cell lysate.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using DNMT1 mouse mAb with DAB staining.

产品名称: DNMT1 Mouse Monoclonal Antibody
产品货号: AMM81255



注意事项 (Note)

For research use only .