

产品名称: THY1 Mouse Monoclonal Antibody
产品货号: AMM81259



产品概述 (Summary)

产品名称 (Production Name)	THY1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,ELISA
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	THY1
别名 (Alternative Names)	CD90
基因 ID (Gene ID)	7070.0
蛋白 ID (SwissProt ID)	P04216.Purified recombinant fragment of human THY1 (AA: 17-132) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	18kDa

研究背景 (Background)

THY1 may play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain <br


产品名称: THY1 Mouse Monoclonal Antibody
产品货号: AMM81259




/>

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using THY1 mouse mAb against T47D (1), HepG2 (2) and PC-12 (3) cell lysate.




Immunofluorescence analysis of HeLa cells using THY1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.


产品名称: THY1 Mouse Monoclonal Antibody
产品货号: AMM81259



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using THY1 mouse mAb with DAB staining.




Immunohistochemical analysis of paraffin-embedded human cerebellum tissues using THY1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rat spleen using THY1 mouse mAb with DAB staining.

产品名称: THY1 Mouse Monoclonal Antibody
产品货号: AMM81259



Immunohistochemical analysis of paraffin-embedded Rat kidney using THY1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .