

产品名称: SYCP3 Mouse Monoclonal Antibody
产品货号: AMM81266



产品概述 (Summary)

产品名称 (Production Name)	SYCP3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SYCP3
别名 (Alternative Names)	COR1; SCP3; SPGF4
基因 ID (Gene ID)	50511.0
蛋白 ID (SwissProt ID)	Q8IZU3. Purified recombinant fragment of human SYCP3 (AA: 27-128) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:200-1:1000, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	27.7kDa

研究背景 (Background)

This gene encodes an essential structural component of the synaptonemal complex. This complex is involved in synapsis,


产品名称: SYCP3 Mouse Monoclonal Antibody
产品货号: AMM81266




recombination and segregation of meiotic chromosomes. Mutations in this gene are associated with azoospermia in males and susceptibility to pregnancy loss in females. Alternate splicing results in multiple transcript variants that encode the same protein.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Immunofluorescence analysis of HepG2 cells using SYCP3 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.




产品名称: SYCP3 Mouse Monoclonal Antibody
产品货号: AMM81266



Flow cytometric analysis of Jurkat cells using SYCP3 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using SYCP3 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human kidney tissues using SYCP3 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .