

产品名称: TFRC Mouse Monoclonal Antibody
产品货号: AMM81272



产品概述 (Summary)

产品名称 (Production Name)	TFRC Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Monkey,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TFRC
别名 (Alternative Names)	T9; TR; TFR; p90; CD71; TFR1; TRFR
基因 ID (Gene ID)	7037.0
蛋白 ID (SwissProt ID)	P02786.Purified recombinant fragment of human TFRC (AA: 608-727) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	85kDa

研究背景 (Background)

Transferrin receptor is a carrier protein for transferrin. It is needed for the import of iron into the cell and is regulated in


产品名称: TFRC Mouse Monoclonal Antibody
产品货号: AMM81272




response to intracellular iron concentration. Low iron concentrations promote increased levels of transferrin receptor, to increase iron intake into the cell. Thus, transferrin receptor maintains cellular iron homeostasis. Expression of human TFR1, but not human TFR2, in hamster cell lines markedly enhanced the infection of viruses pseudotyped with the glycoprotein of Machupo, Guanarito, and Junin viruses, but not with those of Lassa or lymphocytic choriomeningitis viruses. An anti-TFR1 antibody efficiently inhibited the replication of Machupo, Guanarito, Junin, and Sabia viruses, but not that of Lassa virus. TFR1 is a cellular receptor for New World hemorrhagic fever arenaviruses.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using TFRC mouse mAb against Jurkat (1), Hela (2), K562 (3), Cos7 (4), MCF-7 (5), PC-12 (6), NIH/3T3 (7), HEK293 (8), RAJI (9) cell lysate.

注意事项 (Note)

For research use only .