

产品名称: DIS3L2 Mouse Monoclonal Antibody
产品货号: AMM81273



产品概述 (Summary)

产品名称 (Production Name)	DIS3L2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DIS3L2
别名 (Alternative Names)	FAM6A; PRLMNS
基因 ID (Gene ID)	129563.0
蛋白 ID (SwissProt ID)	Q8IYB7 .Purified recombinant fragment of human DIS3L2 (AA: 27-250) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	99.3kDa

研究背景 (Background)

The protein encoded by this gene is similar in sequence to 3'/5' exonucleolytic subunits of the RNA exosome. The exosome


产品名称: DIS3L2 Mouse Monoclonal Antibody
产品货号: AMM81273




is a large multimeric ribonucleotide complex responsible for degrading various RNA substrates. Several transcript variants, some protein-coding and some not, have been found for this gene.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using DIS3L2 mouse mAb against HeLa (1) and HepG2 (2) cell lysate.




Flow cytometric analysis of Jurkat cells using DIS3L2 mouse mAb (green) and negative control (red).

产品名称: DIS3L2 Mouse Monoclonal Antibody
产品货号: AMM81273



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using DIS3L2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using DIS3L2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .