

产品名称: CDX2 Mouse Monoclonal Antibody
产品货号: AMM81289



产品概述 (Summary)

产品名称 (Production Name)	CDX2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDX2
别名 (Alternative Names)	CDX3; CDX-3
基因 ID (Gene ID)	1045.0
蛋白 ID (SwissProt ID)	Q99626.Purified recombinant fragment of human CDX2 (AA: 176-303) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	33.5kDa

研究背景 (Background)

This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded protein is a major


产品名称: CDX2 Mouse Monoclonal Antibody
产品货号: AMM81289




regulator of intestine-specific genes involved in cell growth and differentiation. This protein also plays a role in early embryonic development of the intestinal tract. Aberrant expression of this gene is associated with intestinal inflammation and tumorigenesis.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Hela cells using CDX2 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .