

产品名称: Bpifa2 Mouse Monoclonal Antibody
产品货号: AMM81310



产品概述 (Summary)

产品名称 (Production Name)	Bpifa2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Bpifa2
别名 (Alternative Names)	Psp; Bpifa2e;msplunc2
基因 ID (Gene ID)	19194.0
蛋白 ID (SwissProt ID)	P07743 .Purified recombinant fragment of mouse mSplunc2 (AA: 16-169) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	24.8kDa

研究背景 (Background)


Bpifa2 has strong antibacterial activity against P. aeruginosa.

产品名称: Bpifa2 Mouse Monoclonal Antibody
产品货号: AMM81310




研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Flow cytometric analysis of MCF-7 cells using mSplunc2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using mSplunc2 mouse mAb with DAB staining.

产品名称: Bpifa2 Mouse Monoclonal Antibody
产品货号: AMM81310



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using mSplunc2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .