

产品名称: DAPK3 Mouse Monoclonal Antibody
产品货号: AMM81334



产品概述 (Summary)

产品名称 (Production Name)	DAPK3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DAPK3
别名 (Alternative Names)	ZIP; ZIPK
基因 ID (Gene ID)	1613.0
蛋白 ID (SwissProt ID)	O43293.Purified recombinant fragment of human DAPK3 (AA: 28-161) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	52.5kDa

研究背景 (Background)

Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in

产品名称: DAPK3 Mouse Monoclonal Antibody
产品货号: AMM81334




mammalian cells. These results suggest that DAPK3 may play a role in the induction of apoptosis.


研究领域 (Research Area)

Apoptosis


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using DAPK3 mouse mAb against A431 (1), K562 (2) cell lysate.



Flow cytometric analysis of A431 cells using DAPK3 mouse mAb (green) and negative control (red).

产品名称: DAPK3 Mouse Monoclonal Antibody
产品货号: AMM81334



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using DAPK3 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .