

产品名称: GAB2 Mouse Monoclonal Antibody
产品货号: AMM81338



产品概述 (Summary)

产品名称 (Production Name)	GAB2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GAB2
别名 (Alternative Names)	GAB2
基因 ID (Gene ID)	9846.0
蛋白 ID (SwissProt ID)	Q9UQC2. Purified recombinant fragment of human GAB2 (AA: 332-476) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:200-1:1000, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	74.5kDa

研究背景 (Background)

This gene is a member of the GRB2-associated binding protein (GAB) gene family. These proteins contain pleckstrin


产品名称: GAB2 Mouse Monoclonal Antibody
产品货号: AMM81338




homology (PH) domain, and bind SHP2 tyrosine phosphatase and GRB2 adapter protein. They act as adapters for transmitting various signals in response to stimuli through cytokine and growth factor receptors, and T- and B-cell antigen receptors. The protein encoded by this gene is the principal activator of phosphatidylinositol-3 kinase in response to activation of the high affinity IgE receptor. Two alternatively spliced transcripts encoding different isoforms have been described for this gene.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Immunofluorescence analysis of HeLa cells using GAB2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Secondary antibody from Fisher (Cat#: 35503)


产品名称: GAB2 Mouse Monoclonal Antibody
产品货号: AMM81338



Flow cytometric analysis of Hela cells using GAB2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using GAB2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using GAB2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .