

产品名称: B3GALT4 Mouse Monoclonal Antibody
产品货号: AMM81341



产品概述 (Summary)

产品名称 (Production Name)	B3GALT4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	B3GALT4
别名 (Alternative Names)	GALT2; GALT4; BETA3GALT4
基因 ID (Gene ID)	8705.0
蛋白 ID (SwissProt ID)	O96024.Purified recombinant fragment of human B3GALT4 (AA: 191-359) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	41.5kDa

研究背景 (Background)

This gene is a member of the beta-1,3-galactosyltransferase (beta3GalT) gene family. This family encodes type II


产品名称: B3GALT4 Mouse Monoclonal Antibody
产品货号: AMM81341




membrane-bound glycoproteins with diverse enzymatic functions using different donor substrates (UDP-galactose and UDP-N-acetylglucosamine) and different acceptor sugars (N-acetylglucosamine, galactose, N-acetylgalactosamine). The beta3GalT genes are distantly related to the Drosophila Brainiac gene and have the protein coding sequence contained in a single exon. The beta3GalT proteins also contain conserved sequences not found in the beta4GalT or alpha3GalT proteins. The carbohydrate chains synthesized by these enzymes are designated as type 1, whereas beta4GalT enzymes synthesize type 2 carbohydrate chains. The ratio of type 1:type 2 chains changes during embryogenesis. By sequence similarity, the beta3GalT genes fall into at least two groups: beta3GalT4 and 4 other beta3GalT genes (beta3GalT1-3, beta3GalT5). This gene is oriented telomere to centromere in close proximity to the ribosomal protein S18 gene. The functionality of the encoded protein is limited to ganglioseries glycolipid biosynthesis.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using B3GALT4 mouse mAb against PANC-1 (1), PC-3 (2) cell lysate.

产品名称: B3GALT4 Mouse Monoclonal Antibody
产品货号: AMM81341



Flow cytometric analysis of PANC-1 cells using B3GALT4 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .