

产品名称: TUBE1 Mouse Monoclonal Antibody
产品货号: AMM81362



产品概述 (Summary)

产品名称 (Production Name)	TUBE1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TUBE1
别名 (Alternative Names)	TUBE; dJ142L7.2
基因 ID (Gene ID)	51175.0
蛋白 ID (SwissProt ID)	Q9UJT0.Purified recombinant fragment of human TUBE1 (AA: 314-472) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	53kDa

研究背景 (Background)

This gene encodes a member of the tubulin superfamily. This protein localizes to the centriolar sub-distal appendages that


产品名称: TUBE1 Mouse Monoclonal Antibody
产品货号: AMM81362




are associated with the older of the two centrioles after centrosome duplication. This protein plays a central role in organization of the microtubules during centriole duplication. A pseudogene of this gene is found on chromosome 5.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Flow cytometric analysis of Hela cells using TUBE1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using TUBE1 mouse mAb with DAB staining.

产品名称: TUBE1 Mouse Monoclonal Antibody
产品货号: AMM81362



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using TUBE1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .