

产品名称: MEN1 Mouse Monoclonal Antibody
产品货号: AMM81366



产品概述 (Summary)

产品名称 (Production Name)	MEN1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MEN1
别名 (Alternative Names)	MEAI; SCG2
基因 ID (Gene ID)	4221.0
蛋白 ID (SwissProt ID)	O00255.Purified recombinant fragment of human MEN1 (AA: 392-554) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	68kDa

研究背景 (Background)

This gene encodes menin, a putative tumor suppressor associated with a syndrome known as multiple endocrine neoplasia


产品名称: MEN1 Mouse Monoclonal Antibody
产品货号: AMM81366




type 1. In vitro studies have shown menin is localized to the nucleus, possesses two functional nuclear localization signals, and inhibits transcriptional activation by JunD, however, the function of this protein is not known. Two messages have been detected on northern blots but the larger message has not been characterized. Alternative splicing results in multiple transcript variants.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Hela cells using MEN1 mouse mAb (green) and negative control (red).



产品名称: MEN1 Mouse Monoclonal Antibody
产品货号: AMM81366



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using MEN1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .