

产品名称: RAB27A Mouse Monoclonal Antibody
产品货号: AMM81369



产品概述 (Summary)

产品名称 (Production Name)	RAB27A Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB27A
别名 (Alternative Names)	GS2; RAM; RAB27; HsT18676
基因 ID (Gene ID)	5873.0
蛋白 ID (SwissProt ID)	P51159.Purified recombinant fragment of human RAB27A (AA: FULL(1-221)) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	24.9kDa

研究背景 (Background)

The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. The protein is membrane-bound


产品名称: RAB27A Mouse Monoclonal Antibody
产品货号: AMM81369




and may be involved in protein transport and small GTPase mediated signal transduction. Mutations in this gene are associated with Griscelli syndrome type 2. Alternative splicing occurs at this locus and four transcript variants encoding the same protein have been identified.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissues using RAB27A mouse mAb with DAB staining.



产品名称: RAB27A Mouse Monoclonal Antibody
产品货号: AMM81369



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using RAB27A mouse mAb with DAB staining.

注意事项 (Note)

For research use only .