

产品名称: RHOF Mouse Monoclonal Antibody
产品货号: AMM81377



产品概述 (Summary)

产品名称 (Production Name)	RHOF Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RHOF
别名 (Alternative Names)	RIF; ARHF
基因 ID (Gene ID)	54509.0
蛋白 ID (SwissProt ID)	Q9HBH0.Purified recombinant fragment of human RHOF (AA: 1-84) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	23.6kDa

研究背景 (Background)

RHOF is a protein-coding gene. Diseases associated with RHOF include mycosis fungoides, and lung cancer, and among its


产品名称: RHOF Mouse Monoclonal Antibody
产品货号: AMM81377




related super-pathways are Signaling by Rho GTPases and Signaling by GPCR. GO annotations related to this gene include GTP binding and GTPase activity. An important paralog of this gene is RHOB.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HeLa cells using RHOF mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .