

产品名称: VTN Mouse Monoclonal Antibody
产品货号: AMM81386



产品概述 (Summary)

产品名称 (Production Name)	VTN Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VTN
别名 (Alternative Names)	VN; V75; VNT
基因 ID (Gene ID)	7448.0
蛋白 ID (SwissProt ID)	P04004.Purified recombinant fragment of human VTN (AA: 20-199) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	54.3kDa

研究背景 (Background)

The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell


产品名称: VTN Mouse Monoclonal Antibody
产品货号: AMM81386




adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Hela cells using VTN mouse mAb (green) and negative control (red).




Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using VTN mouse mAb with DAB

产品名称: VTN Mouse Monoclonal Antibody
产品货号: AMM81386



staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using VTN mouse mAb with DAB staining.

注意事项 (Note)

For research use only .