

产品名称: APP Mouse Monoclonal Antibody
产品货号: AMM81396



产品概述 (Summary)

产品名称 (Production Name)	APP Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APP
别名 (Alternative Names)	AAA; AD1; PN2; ABPP; APPI; CVAP; ABETA; PN-II; CTFgamma
基因 ID (Gene ID)	351.0
蛋白 ID (SwissProt ID)	P05067.Purified recombinant fragment of human APP (AA: 483-699) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	87kDa

研究背景 (Background)

This gene encodes a cell surface receptor and transmembrane precursor protein that is cleaved by secretases to form a

产品名称: APP Mouse Monoclonal Antibody
产品货号: AMM81396




number of peptides. Some of these peptides are secreted and can bind to the acetyltransferase complex APBB1/TIP60 to promote transcriptional activation, while others form the protein basis of the amyloid plaques found in the brains of patients with Alzheimer disease. Mutations in this gene have been implicated in autosomal dominant Alzheimer disease and cerebroarterial amyloidosis (cerebral amyloid angiopathy). Multiple transcript variants encoding several different isoforms have been found for this gene.


研究领域 (Research Area)

Apoptosis, Notch signaling pathway

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Hela cells using APP mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .