

产品名称: APBB1IP Mouse Monoclonal Antibody
产品货号: AMM81406



产品概述 (Summary)

产品名称 (Production Name)	APBB1IP Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APBB1IP
别名 (Alternative Names)	RIAM; INAG1; PREL1; RARP1
基因 ID (Gene ID)	54518.0
蛋白 ID (SwissProt ID)	Q7Z5R6 .Purified recombinant fragment of human APBB1IP (AA: 1-151) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	73.2kDa

研究背景 (Background)

APBB1IP (amyloid beta (A4) precursor protein-binding, family B, member 1 interacting protein) is a protein-coding gene.


产品名称: APBB1IP Mouse Monoclonal Antibody
产品货号: AMM81406




Diseases associated with APBB1IP include alzheimer's disease, and melanoma, and among its related super-pathways are p130Cas linkage to MAPK signaling for integrins and Platelet Aggregation (Plug Formation). GO annotations related to this gene include phospholipid binding. An important paralog of this gene is GRB7.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Hela cells using APBB1IP mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human lymph tissues using APBB1IP mouse mAb with DAB staining.

产品名称: APBB1IP Mouse Monoclonal Antibody
产品货号: AMM81406



注意事项 (Note)

For research use only .