

产品名称: TFF2 Mouse Monoclonal Antibody
产品货号: AMM81418



产品概述 (Summary)

产品名称 (Production Name)	TFF2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TFF2
别名 (Alternative Names)	SP; SML1
基因 ID (Gene ID)	7032.0
蛋白 ID (SwissProt ID)	Q03403.Purified recombinant fragment of human TFF2 (AA: 20-125) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	14.3kDa

研究背景 (Background)

Members of the trefoil family are characterized by having at least one copy of the trefoil motif, a 40-amino acid domain


产品名称: TFF2 Mouse Monoclonal Antibody
产品货号: AMM81418




that contains three conserved disulfides. They are stable secretory proteins expressed in gastrointestinal mucosa. Their functions are not defined, but they may protect the mucosa from insults, stabilize the mucus layer and affect healing of the epithelium. The encoded protein inhibits gastric acid secretion. This gene and two other related trefoil family member genes are found in a cluster on chromosome 21.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Hela cells using TFF2 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .