

产品名称: HH3 Mouse Monoclonal Antibody
产品货号: AMM81425



产品概述 (Summary)

产品名称 (Production Name)	HH3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HH3
别名 (Alternative Names)	HIST3H3;H3t; H3.4; H3/g; H3FT
基因 ID (Gene ID)	8290.0
蛋白 ID (SwissProt ID)	Q16695.Synthesized peptide of human HH3 (AA: ARTKQTAR(AcK)STG-C).

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	15.5kDa

研究背景 (Background)

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs


产品名称: HH3 Mouse Monoclonal Antibody
产品货号: AMM81425




of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a member of the histone H3 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is located separately from the other H3 genes that are in the histone gene cluster on chromosome 6p22-p21.3.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using AHH3 mouse mAb against NIH3T3 (1), HeLa (2), K562 (3) cell lysate.

产品名称: HH3 Mouse Monoclonal Antibody
产品货号: AMM81425



Flow cytometric analysis of NIH/3T3 cells using HH3 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .