

产品名称: WTAP Mouse Monoclonal Antibody
产品货号: AMM81445



产品概述 (Summary)

产品名称 (Production Name)	WTAP Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	WTAP
别名 (Alternative Names)	PNAS-132;hFL(2)D;KIAA0105
基因 ID (Gene ID)	9589.0
蛋白 ID (SwissProt ID)	Q15007.Purified recombinant fragment of human WTAP (AA: 91-201) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	44.2kDa

研究背景 (Background)

The Wilms tumor suppressor gene WT1 appears to play a role in both transcriptional and posttranscriptional regulation of


产品名称: WTAP Mouse Monoclonal Antibody
产品货号: AMM81445




certain cellular genes. This gene encodes a WT1-associated protein, which is a ubiquitously expressed nuclear protein. Like WT1 protein, this protein is localized throughout the nucleoplasm as well as in speckles and partially colocalizes with splicing factors. Alternative splicing of this gene results in several transcript variants encoding three different isoforms.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Flow cytometric analysis of Hela cells using WTAP mouse mAb (green) and negative control (red).



产品名称: WTAP Mouse Monoclonal Antibody
产品货号: AMM81445



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using WTAP mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using WTAP mouse mAb with DAB staining.

注意事项 (Note)

For research use only .